

Flow cytometry Lab. Royan Institute for Stem cell Biology and Technology



Direct Intracellular Staining



5-10x105 cells / test

Detach cells from culture dish or flask

Wash with PBS⁻(1500 rpm for 5min) at 4°C, discard supernatant

Add 100 μ l cold paraformaldehyde 2-4%; incubate for 20 min at 4°C

Wash cells with PBS⁻/Tween-0/05% (1500 rpm for5 min), discard supernatant

Permeabilization by adding 100 µl (0.1-0.2%) Triton X -100, incubate for 10min at RT

Wash cells with PBS⁻/Tween-0/05% (1500 rpm for5 min), discard supernatant

Add 500 μl secondary host serum 10% (in lack of host serum, use BSA 5%) mix and incubate for 20 min at RT or 30 min at 4°C

Wash cells with PBS /Tween-0/05% (1500 rpm for5 min)

Isolate cells pellet for Unstain control, Isotype control and tests in separated tubes (each tube consist of 100 μ l perm/wash buffer)

Add fluorescent conjugated Antibody directly (datasheet's recommendation), mix and incubate for 30-45 min at 4°C in dark Add isotype control Antibody (datasheet's recommendation), mix and incubate for 30-45 min at 4°C in dark

Rinse the cells with PBS-/Tween-0/05% (1500 rpm for5 min) 2 times, remove the supernatant, and resuspend cells with 100 μl PBS⁻/Tween-0/05%

Run and analysis on Flow Cytometer or fix samples by adding 100ul paraformaldehyde 1% and keep at 4°C in dark up to 24h.